

316L SST Optional 100C UL Certified 16 CF Vacuum Oven With 6 Heating Shelves

\$13,421.50 \$15,790.00



<u>Description</u> <u>Specifications</u> <u>Standard Packages</u>

Description

Specifications

Vacuum oven specifications

Electrical requirements	110V +/-10% 60Hz 1-PH 13.6A, 1,500 watts or 220V +/-10% 50/60Hz 1-PH 6.8A, 1,500 watts
Heating	250 watts x 6, built into 6 aluminum shelves, for the fastest, most energy efficient and uniform heating.
Chamber	Dimensions: 25.6 x 30 x 35.5" (WxDxH) Material: Stainless steel Capacity: 16 cubic feet
Temperature control We use cookies to make your experience better. To comply with the the cookies. Learn more. Allow Cookies	Controller: low proportional gain, microcomputer PID controlled with LCD display e new e-Privacy directive we need to ask for your consent to set Range: ambient to 100°C (or 212°F) Display units: Fahrenheit or Celsius Uniformity: +/- 5% of setpoint Dwelling timer range: 1 to 9999 minutes

	Mechanical vacuum gauge range: 0 to 30 inch mercury Mechanical vacuum gauge type: oil-filled Vacuum port: KF25 flange x 1 Vent port: KF25 flange x 1 Vacuum pump: Sold separately Approximate time to reach full vacuum: With Edwards E2M28 21 cfm pump: 8 minutes With Ai SuperVac 11.3 cfm pump: 11 minutes
Weight	Unit: 710 Lbs, shipping: 1114 Lbs
Dimensions (WxDxH)	Unit: 33 x 37 x 52" Shipping: 46 x 39 x 60"
Shelves	Heated shelves: 6 non-removable aluminum shelves, 25 x 29.25 inches (WxD) each Optional stackable aluminum shelves: 25 x 29.25 x 2.375 inches (WxDxH) Distance between shelves (without stackable shelves): 5 inches Distance between shelves (with 6 stackable shelves): 2.375 inches Total area without stackable shelves: 731 x 6 = 30 sq ft or 4386 sq inch Total area with 6 stackable shelves: 731 x 12 = 60 sq ft or 8772 sq inch Max. load capacity per shelf: 20 Lbs
Door gasket material	Silicone or Viton
Observation window	1/2" tempered safety glass
Safety	Built-in circuit breaker Overheat shutoff protection Secondary over-temp protection dial
In-door lights	White LED
Compliance We use cookies to make your experience better. To comply with the war cookies. Learn more.	UL (E482564), CE e new e-Privacy directive, we need to ask for your consent to set 2 years

Product Video

Across International AT Series Vacuum Ovens Key features



AT vacuum oven installation and calibration (English)

Want to know how to install the Across International Vacuum Oven in your Lab?



AT vacuum oven installation and calibration (Spanish)

Ai AccuTemp 1.9 Cu Ft Vacuum Oven Installation Calibration (Spanish Version)



We use cookies to make your experience better. To comply with the new e-Privacy directive, we need to ask for your consent to set the cookies. <u>Learn more</u>.