Semi-Automatic Production Scale Capsule Filling Machine





The cost-effective, low maintenance leader

CAP8[®] Capsule Filling Machine

Our easy to operate CAP8^{*} is a low-maintenance and cost-effective filling machine that offers semi-automatic operation and reliable performance. The CAP8 is designed to allow you to scale up to production-level volumes without the need for substantial investment.

Specifications

Production Capacity^{*}

Speed per hour
17,000–25,000
20,000–29,000
17,000–25,000
14,000–22,000
9,000–15,000
7,000–11,000

AAA**, AA**, AAel**, A, B and C: Contact Capsugel's Technical Support Group for details.

* Rates depend on operator's skill level and density and type of material being filled. ** Conversion kit required.





Dimensions 136 cm x 77 cm x 191 cm high

Weight 440 kg (970 lbs)

Power Requirements 200-240 VAC 3-Phase on a 20-amp circuit; WYO or Delta Power

Other



All orders include 1 set of change parts and choice of auger (Spade, One-Way, Two-Way, or Three-Way).

To bring your production on line immediately, we offer an installation package including two days of technical training and up to three days of operator training.

User-friendly Features

A wide array of innovative features ensures years of efficient filling and reliability.

- Table-top pneumatic closing station for easy integration of capsule finishing equipment
- Powder hopper is removable for cleaning
- Easy auger adjustment to change fill weights
- Easy closed capsule length adjustment
- Quick and accurate change part set-up
- "Toolless" changeover of rectifier parts
- Easy speed selection
- Large powder recovery station
- Robust electronic controls
- Easy to remove rotary table and ring carrier assemblies

GMP Features

The CAP8's advanced technology and components ensure cleanliness and protect product integrity.

- All contact surfaces are made with FDA approved materials
- Stainless steel ring carrier, rotary table and powder hopper
- The internal components of the cabinet are protected against powder infiltration
- Easy to clean

8

Safety Features

The unique safety features of the CAP8 protect operators and mechanics.

- Remote E-stop for operator safety
- Guards and proximity switch on closing station
- Rectifier guard interlocks
- CE compliant

Optional Features

- Casters
- Fill ring rest tray
- Free-standing vibrator

Support

Capsugel's highly experienced technical service team is always ready to assist customers. Offering hands-on support for shop-floor operators and managers, they can optimize production performance and help prevent problems before they arise. They will also work closely with your key production and maintenance people to identify areas of improvement and to develop machine-specific action plans.



For more information on the CAP8® Capsule Filling Machine and how it can help improve your efficiency, contact your Capsugel sales representative (see list of worldwide locations on back) or visit us on the web at www.capsugel.com.

💷 n 🛄

00

AMERICAS

Greenwood, South Carolina Capsugel 535 North Emerald Road Greenwood, SC 29646 USA Tel: 888-783-6361 or 864-223-2270 Fax: 888-783-6360

Sales Office México City, México

Pfizer Consumer Health Care México, S. de R.L. de CV. Div. Capsugel Paseo de los Tamarindos #40 Col. Bosques de las Lomas 05120 México, D.F. México Tel: (5255) 9177 3120 Fax: (5255) 5257 5174

Sales Office Rio de Janeiro, Brazil

Laboratórios Pfizer Ltda. Div. Capsugel Avenida Das Américas #8445 Salas 1202 to 1207 Barra da Tijuca CEP: 22793-081 Rio de Janeiro, RJ. Brazil Tel: +55-21-2432 5266 Fax: +55-21-2432-5270

EUROPE

Bornem, Belgium Capsugel Bornem Rijksweg 11 B-2880 Bornem Belgium Tel: 32-(0)3-889-05-11 Fax: 32-(0)3-889-26-22

Colmar, France Capsugel Colmar 10, Rue Timken F-68027 Colmar, Cedex France Tel: 33-(0)3-89-20-57-09 Fax: 33-(0)3-89-41-48-11

 Ploërmel, France

 Capsugel

 Z.I. Camagnon - BP 320

 56 803 Ploërmel, Cedex

 France

 Tel: 33-(0)2-97-72-10-00

 Fax: 33-(0)2-97-72-15-00

Cambridge, UK Capsugel Cambridge Research Centre Compass House, Vision Park Chivers Way, Histon, Cambridge CB4 9ZR United Kingdom Tel: 44-(0)-1223-598480 Fax: 44-(0)-1223-598070

For information on our other machines and services designed to speed your pharmaceutical development time, such as the Xcelodose[™]Precision Micro-filling System, the CFS 1000 and the LEMS 60[™], please visit us at www.capsugel.com.



Ayutthaya, Thailand Capsugel (Thailand) Co., Ltd. 1/82 Moo 5, Rojana Industrial Park Rojana Road, Tambol Khanham, Amphur U-Thai Ayutthaya 13210, Thailand Tel: 66 35 33 4000 Fax: 66 35 22 7088

Jakarta, Indonesia

ASIA

P.T. Capsugel Indonesia JI. Raya Bogor Km 42 P.O. Box 15/CBI Cibinong 16916, Jawa Barat Indonesia Tel: 62-21-875-2226 Fax: 62-21-875-2227

Suzhou, China

Suzhou Capsugel Ltd. No. 369 Suhong Middle Road, Suzhou Industrial Park Jiangsu Province, P.R.China, 215021 Tel: 86-512-6258-5188 Fax: 86-512-6258-5188

Sales Office Hong Kong

16/F, Stanhope House, 738 Kings Road North Point, Hong Kong Tel: (852) 2811 9711, 2562 5567 (direct) Fax: (852) 2214 1737

Sagamihara, Japan

Capsugel Japan KK 3-36 Minamihashimoto, 4-chome Sagamihara City Kanagawa Prefecture 229-1133 Japan Tel: 81-427-00-6700 Fax: 81-427-00-6719

Sydney, Australia

Capsugel 38-42 Wharf Road West Ryde NSW 2114 Australia Tel: 61-2-9813-4405 Fax: 61-2-9813-4411

