








Holder Type	Description	Applications	Maximum Size mm (in) WxHxD	Exposure Size mm (in) WxH	Instrument	Rack Capacity
 SL-3T Single exposure window with spring clip back Part Number 19163900	General purpose holder for exposure of thin specimens	Textiles, plastic film, automotive interior	67 x 145 x 3 (2.6 x 5.7 x 0.1)	50 x 121 (2.0 x 4.7)	Ci3000+ Ci4000 Ci5000	20 68 111
 SL-3T with Glass Single exposure window with glass and adjustable back Part Number 07303900	Exposure through secondary filter to closely simulate end-use environment	Textiles, paper, plastic film, carpet, automotive interior	67 x 145 x 15 (2.6 x 5.7 x 0.6)	50 x 121 (2.0 x 4.7)	Ci3000+ Ci4000 Ci5000	20 68 111
 CD-3T Three exposure windows with spring clip back Part Number 20215700	General purpose holder with three exposure windows for thin specimens	Textiles, paper, plastic film, automotive interior	67 x 145 x 3 (2.6 x 5.7 x 0.1)	3 windows: 38 x 50 (1.5 x 2.0)	Ci3000+ Ci4000 Ci5000	20 68 111
 CD-3T with Glass Three exposure windows with glass, spring clip back Part Number 07303800	Exposure through secondary filter to closely simulate end-use environment	Textiles, paper, plastic film, wood, automotive interior	67 x 145 x 15 (2.6 x 5.7 x 0.6)	3 windows: 38 x 50 (1.5 x 2.0)	Ci3000+ Ci4000 Ci5000	20 68 111
 RD-3T Single window or 3-window exposure with "bulldog" clip holder Part Number 20017900	Versatile holder to accommodate a range of specimen widths and thicknesses	Coatings on various substrates, plastics, textiles, glass	77 x 152 x 10 (3.0 x 6.0 x 0.4)	57 x 134 (2.3 x 5.3)	Ci3000+ Ci4000 Ci5000	20 68 111
 TEX-3T with Mask Single exposure window with mask, adjustable Part Number 19186700	Holder with masked exposure and adjustable back to allow for thick sample	Textiles, foam, foam-backed materials	45 x 134 x 12 (1.7 x 5.3 x 0.5)	19 x 119 (0.7 x 4.7)	Ci3000+ Ci4000 Ci5000	29 104 170
 Panel Holder Part Number 19188501	Holds 76 x 152 (3 x 6) panels for exposure tests of coatings	Coatings, rigid plastic, wood	76 x 152 x 9 (3.0 x 6.0 x 0.37)	76 x 146 (3.0 x 5.7)	Ci3000+ Ci4000 Ci5000	16 56 87
 Drop-in Specimen Bar Holder Part Number 19184600	Holds vertical specimen bars	Plastics	77 x 144 x 3 (3.0 x 5.6 x 0.1)	76 x 125 (3.0 x 4.9)	Ci3000+ Ci4000 Ci5000	15 50 84

This chart is a representative sample of specimen holders available for a variety of Atlas® instruments. For specific information about specimen holders that best meet your needs, please contact your local Atlas representative.



Specimen Holders

Accelerated Weathering Products

Rack Capacity	Instrument	Exposure Size mm (in) WxH	Maximum Size mm (in) WxHxD	Applications	Description	Holder Type
11	Xenotest® Alpha+ Xenotest 150 S+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Textiles, plastics, coatings, paper	Holds specimen up to 3 mm thick	Regular Specimen Holder Part Number 56075142
11	Xenotest Alpha+ Xenotest 150 S+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Carpets, foam, foam-backed materials, thick panels	Holds specimen up to 15 mm thick	Specimen Holder for Thick Specimens Part Number 56077987
only 1 needed (replaces one holder)	Xenotest Alpha+ Xenotest 150 S+	—	135 x 45 (5.3 x 1.8)	Blue Scale fabric	Holder for Blue Scale fabric during weathering tests	Specimen Holder for Blue Scale Part Number 56050873
16	Xenotest Beta+	295 x 70 (11.6 x 2.8)	310 x 80 (12.2 x 3.2)	Textiles, plastics, coatings, paper, window profiles	Holds specimen up to 10 mm thick	Standard Specimen Holder (including back plate) Part Number 56076543
16	Xenotest Beta+	90 x 58 (3.5 x 2.3)	100 x 68 (3.9 x 2.7)	Textiles, plastics, coatings, paper	Holds 3 specimens up to 10 mm thick, fixed by a back plate and a cover mask	Special Specimen Holder, 3 Segments Part Number 56076699
28 19	Xenotest Beta+ Xenotest 220/220+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Carpet, foam-backed textiles	Holds 2 specimens up to 10 mm thick, fixed by 2 cover masks	Special Specimen Holder, 2 Segments Part Number 56076700
28 19	Xenotest Beta+ Xenotest 220/220+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Textiles, paper, leather	Holds 2 samples up to 3 mm thick, fixed by 2 cover masks	Specimen Holder Standard, 2 Segments Quick-insertion Part Number 56079776
28 19	Xenotest Beta+ Xenotest 220/220+	270 x 45 (10.6 x 1.8)	320 x 48 (12.6 x 1.9)	Technical textiles	Holds 1 specimen up to 3 mm thick	Specimen Holder Special TS, 1 Segment Quick-insertion, Part Number 56079797
only 1 needed (replaces one holder)	Xenotest 220/220+	—	135 x 45 (5.3 x 1.8)	Blue Scale fabric	Holds Blue Scale fabric for weathering tests	Specimen Holder for Blue Scale (Special holder kit, 2 Segments) Part Number 56078036
15	SUNTEST® XLS+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Plastics, coatings, paper, leather	General purpose holder for various materials	Set of 15 Specimen Trays Part Number 56079636
34	SUNTEST XXL/XXL+	121 x 35 (4.8 x 1.4)	135 x 45 (5.3 x 1.8)	Plastics, coatings, paper	General purpose holder for various materials	Set of Five Specimen Trays for 7 Specimens Part Number 56079259
8	SUNTEST XXL/XXL+	295 x 70 (11.6 x 2.8)	310 x 80 (12.2 x 3.2)	Plastics, coatings, paper	General purpose holder for various materials	Set of Two Specimen Trays for 4 Specimens Part Number 56079261

This chart is a representative sample of specimen holders available for a variety of Atlas® instruments. For specific information about specimen holders that best meet your needs, please contact your local Atlas representative.