


# MediSeal

Pharma Packaging Technology


### MediSeal – synergies under one strong roof

MediSeal belongs to the KÖRBER **MEDIPAK** group of companies and – together with the Swiss packaging specialists Rondo and Dividella – offer a full range of packaging systems. Exploit the synergies of a strong group and benefit from low-cost packaging solutions for your production!



# MediSeal

Pharma Packaging Technology



Small and perfectly formed.  
The CP2 blister machine

CP2  
Technical Data

## CP2 blister machine

max. mm

### Advantages:

- Small machine size (from 3 m)
- Compact tools
- Simplicity of operation
- No limitations on choice of size and foil

Foil width

Draw-off length

Useable area

Form depth

Roll diameter

Base foil

Cover foil

Output, blisters/minute

150

90 (100\*)

142 x 88 (142 x 98\*)

12 (17\*)

400

220

120

\*special version

Subject to technical modifications

### Other KÖRBER MEDIPAK companies

#### MediSeal GmbH

Flurstraße 65  
D-33758 Schloss Holte  
Germany  
T +49 (0)52 07 888-0  
F +49 (0)52 07 888-299  
info@mediseal.de  
www.mediseal.de

#### Dividella AG

Werdenstrasse 76  
CH-9472 Grabs  
Switzerland  
T +41 (0)81 750 33 66  
F +41 (0)81 750 33 43  
info@dividella.ch  
www.dividella.com

#### Rondo AG

Gewerbestrasse 11  
CH-4123 Allschwil  
Switzerland  
T +41 (0)61 486 87 87  
F +41 (0)61 486 87 50  
sales@rondodruck.ch  
www.rondodruck.ch

#### KÖRBER MEDIPAK sales and service companies

**KÖRBER MEDIPAK NA INC.** 14501, 58th Street North • Clearwater, FL 33760 • USA  
T +1 727 5 32 65 10 • F +1 727 5 32 65 37 • info@kmedipak.com

**KÖRBER MEDIPAK France Sarl** Energy Park 4 • 155-159 rue du Dr. Bauer • 93400 Saint Ouen • France  
T +33 1 40 11 40 11 • F +33 1 40 11 51 10 • info@fr.koerber-medipak.com


**KÖRBER MEDIPAK UK & RoI** c/o MediSeal GmbH • Flurstraße 65 • D-33758 Schloß Holte • Germany  
T +44 (7841) 678 122 • info@uk.koerber-medipak.com

**KÖRBER MEDIPAK Asia-Pacific** c/o Hauni Singapore Pte. Ltd. • 6 Temasek Boulevard  
#30-04 Suntec City Tower 4 • Singapore 038986  
T +65 68 35 98 86 • F +65 68 36 92 96 • info@koerber-medipak.com

www.koerber-medipak.com



MediSeal has acquired outstanding know-how related to processing aluminum base foils.





## CP2 blister machine

3 metres in length,  
up to 120 blisters per minute.

The entire upper part of the station can be raised for simple tool changes in the forming station.

The CP2 can be fitted out with a variety of dosing solutions to accommodate all products.



CP2

# The MediSeal CP2

## Flexibility in the smallest possible space.

*The CP2 is a very compact machine of table construction with minimum space requirements and an output of up to 120 blisters per minute.*

*The machine, which has been on the market for more than 15 years, has proved itself in the packaging of small lot sizes, special applications and samples for clinical trials.*

### 3 metres for your production.

The CP2 is equipped as standard with platen sealing. This allows processing of many different foil types, such as paper-coated aluminum foil for childproof packs in the cover foil area as well as all known base foils such as PVC, PVDC, aluminum PP or Aclar.

Thanks to the flexible pocket design, it is possible to package special product such as medical devices, ampoules or cosmetic products with no problems.

Dosing can be adapted to many different production requirements. For example, automatic dosing solutions are available for all known products. Naturally, manual insertion is also possible – if necessary from both sides.

The range of sizes allows the most compact tool sizes, hence the low space requirements for the tooling and the simplest handling of parts during size changes.

The comparatively high output of the CP2 is also suitable for the constantly increasing quantities needed for clinical trial samples.